# RETROGRADE IR PLATFORM

AM0400 - AM0410 - AM0420

## INSTRUCTIONS FOR USE

(Instrucciones de uso / Mode d'emploi / Istruzioni per l'uso / Gebrauchsanweisung / Gebruiksaanwijzing)

Issue Date: 05/2021

Revision: A

EC REP

### EC Representative

Emergo Europe  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

### Manufacturer

Adept Medical Ltd  
2-6 McDonald St,  
Morningside  
Auckland 1025, New Zealand

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

	<b>A</b>	Retrograde Main Platform x1
	<b>B</b>	Front Base x1
	<b>C</b>	Rear Base x1
	<b>D</b>	Leg x 4
	<b>E</b>	Leg Extension x 2
	<b>F</b>	Extension* x1

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















# Setup

(Instalación / Installation / Configurazione / Einrichtung / Installatie)

## 5 Unlock the leg to adjust the platform height

- / Desbloquee la pata para ajustar la altura de la plataforma
- / Déverrouillez le pied pour régler la hauteur de la plateforme
- / Sbloccare la gamba (Leg) per regolare l'altezza della piattaforma
- / Entriegeln Sie den Fuß zum Einstellen der Plattformhöhe
- / Ontgrendel de poot (leg) om het platform in hoogte te verstellen

## 6 ≤ 10 kg

## Spare Parts

(Piezas / Pièces de rechange / Parti di ricambio / Ersatzteile / Reserveonderdelen)

### Product Code

(Código de producto / Code produit / Codice prodotto / Produktcode / Productcode)

	Retrograde IR Platform without Extension	AM0400
	IR Platform Extension	AM0410
	Retrograde IR Platform + Extension	AM0420

### Component Codes

(Código de componentes / Codes des composants / Codici dei componenti / Komponentencodes / Onderdeelcodes)

	<b>A</b>	Retrograde Main Platform	M3660
	<b>B</b>	Front Base	M3662
	<b>C</b>	Rear Base	M3663
	<b>D</b>	Leg	M3664
	<b>E</b>	Leg Extension	M3693
	<b>F</b>	Extension	AM0410

## Image Table Compatibility

(Compatibilidad de la mesa de imagen / Compatibilité de la table d'examen / Compatibilità del tavolo per imaging / Kompatibilität Bildgebungstische / Compatibiliteit van MRI-tafel)

### Ensure table geometry falls within this specification

- / Asegúrese de que la geometría de la mesa cumpla estas especificaciones
- / Assurez-vous que la géométrie de la table correspond à cette spécification
- / Accertarsi che la geometria del tavolo rientri in questa specifica
- / Stellen Sie sicher, dass die Geometrie des Bildgebungstischs innerhalb dieser Angaben liegt
- / De geometrie van de tafel moet binnen deze specificatie vallen

<b>Table Width:</b> / Anchura de la mesa / Largeur de la table / Larghezza tavolo / Tischbreite / Tafelbreedte	45 cm - 70 cm
<b>Table Length:</b> / Longitud de la mesa / Longueur de la table / Lunghezza tavolo / Tischlänge / Tafellengte	≥240 cm
<b>Table Angulation:</b> / Angulación de la mesa / Angulation de la table / Inclinazione tavolo / Winkelstellung des Tisches / Tafelhoek	+/- 15° lateral & longitudinal +/- 15° lateral y longitudinal +/- 15° latéral et longitudinal +/- 15° laterale e longitudinale +/- 15° lateral und longitudinal +/- 15° lateraal en longitudinaal
<b>Table Surface:</b> / Superficie de la mesa / Surface de la table / Superficie del tavolo / Tischoberfläche / Tafeloppervlak	Flat Plana Plane Piatto Flach Vlak

# ESSENTIAL INFORMATION

Información básica / Informations essentielles / Informazioni essenziali / Grundlegende Informationen / Essentiële informatie

## 11 Essential Information (EN)

- Information
- Initial Checks & Cleaning
- Disinfection
- Disinfection Instructions
- Approved Cleaners List
- Disposal

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## Information

### Intended Purpose

To provide a work surface over a supine patient for femoral transarterial diagnosis and intervention via Retrograde approach on an Angiography imaging table.

### Intended User/Training Requirement

Intended to be used by trained medical professionals.

### Patient Target Group

Adult male or female patient requiring diagnostic imaging and image guided therapy.

### Contraindications

- Not to be used with patients over 135Kg.
- Not for paediatric use.

### Warnings and Cautions

- Ensure IFU is read prior to use.
- Ensure Preparatory Cleaning is conducted.
- Ensure Preparatory Product Check is conducted.
- Ensure all items placed on the Retrograde Main Platform do not exceed a total weight of 10Kg.
- Do not use with patients over 135Kg.
- Not for paediatric use.

### Incident Reporting Guidelines

For product complaints and incidents, please complete form on

[adeptmedical.co.nz/repairs](https://adeptmedical.co.nz/repairs)

## Initial Checks & Cleaning

### Cleaning

Refer to Disinfection Instructions.

### Product Checks

- No exposed Carbon Fibre.
- No visible damage or sharp edges, e.g. cracks.
- **Leg:** rotate all four blue locking collars to ensure they are free to rotate and smoothly lock and unlock.
- **Leg:** ensure that the Leg is not free to move up and down in the locked position.
- **Leg:** ensure that the Legs can move up and down when pushed in the unlocked position.
- **Leg:** Ensure the Legs are securely assembled on the Retrograde Main Platform (refer to page 5).
- **Leg:** Ensure Legs are secured in the Lower Bases (refer to page 5).
- **Leg:** Ensure Leg is locked during procedure (refer to page 7).
- **Extension:** Ensure the Extension is fully inserted on the Retrograde Main Platform and the Extension Clips are fully locked (refer to page 6).

## Disinfection

### Warnings

- Insufficient cleaning may compromise disinfection process and lead to transmission of pathogens.
- Do not use any disinfectant products or abrasive/corrosive agents which are not on the approved cleaners list.
- Always read manufacturer's instructions and consult the manufacturers MSDS for cleaning and disinfectant products.
- Do not use any cleaning procedure which exceeds 65°C/149°F.
- Some disinfectants may cause slight discolouration to the blue material used on some components within the product range. This will not affect the strength and the product will remain fit for purpose.

### Limitations on Processing

Discontinue use if:

- Any cracks or breakages are present.
- Any exposed carbon fibre is visible.

## Disinfection Instructions

### Initial treatment at the point of use

It is important to clean the product once it is removed from the packaging and after each procedure. Ensure all areas including joints, clips, sockets, brackets and levers are thoroughly cleaned after use to remove all contaminant build-up that may be present following a procedure.

### Preparation before cleaning

- Remove the product from the table.
- Disassemble Retrograde IR Platform into Leg, Leg Extension, Rear Base, Front Base, Extension (where applicable) and Retrograde Main Platform.
- Ensure the Legs are fully extended.

[adeptmedicaltraining.com/productcleaning/Retrograde-IR-Platform](https://adeptmedicaltraining.com/productcleaning/Retrograde-IR-Platform)

### Cleaning: Manual

- If any visible residue is present, rinse under water, but avoid submerging the device.
- A soft bristle brush can be used.
- Ensure the device is completely dried prior to disinfection.

### Disinfection

Refer to the Approved Cleaners List (following page) for selection of an appropriate disinfecting agent.

#### Retrograde Main Platform:

- With a cloth dampened by an approved disinfectant, wipe all surfaces, including the threaded holes on the underside.

#### Extension (when applicable)

- With a cloth dampened by an approved disinfectant, wipe all surfaces, including entirety of slot and around the blue locking clips.

#### Front and Rear Base:

- With a cloth dampened by an approved disinfectant, wipe all surfaces, including leg sockets and blue locking clips.

#### Leg:

- With a cloth dampened by an approved disinfectant, wipe all surfaces.

#### Leg Extension:

- With a cloth dampened by an approved disinfectant, wipe all surfaces, including entirety of the thread.

### Drying

All components should be dried thoroughly before use.

## Approved Cleaners List

### Approved Disinfectant by Brands

#### Liquids:

ORION Laboratories  
70% Isopropyl Alcohol

Jaychem Industries  
2% Chlorhexidine / 70% Alcohol

Betadine  
Povidone Iodine (7.5%)

BODE Chemie GmbH  
Kohrsolin FF

STERIS  
Coverage Spray TB

CaviCide  
Metrex Research

#### Dissolvable:

Du Pont  
Rely+On Virkon Tablets

#### Surface Wipes:

PDI Sani-Cloth Bleach

PDI Sani-Cloth AF3

Reynard Health Supplies  
Surface Disinfectant Wipes

Metrex Research  
CaviWipes

Clinell Wipes

Clorox Hydrogen Peroxide Wipes

### Approved Active Disinfectant Ingredients

#### Acid Based:

≤10% Malic Acid CAS 6915-15-7  
≤6% Sulfamic acid CAS 5329-14-6

#### Alcohol Based:

≤5% 2-Butoxyethanol CAS 111-76-2  
≤10% Butyldiglycol CAS 112-34-5  
≤70% Isopropyl Alcohol (Propan-2-ol) CAS 67-63-0  
≤10% Tridecanol CAS 69011-36-5  
<10% Alcohols, C12-14, ethoxylated CAS 68439-50-9

#### Sulphate Based:

≤0.1% PHMB CAS 27083-27-8  
≤55% Potassium Peroxymonosulfate CAS 70693-62-8  
≤3% Potassium Persulfate CAS 7727-21-1

#### Ammonium and Chloride Based:

≤5% Benzalkonium Chloride CAS 68424-85-1  
≤0.28% Benzethonium Chloride CAS 121-54-0  
≤10% Benzyl-C23-18-Alkyl-dimethyl Ammonium Chloride CAS 8001-54-5  
≤2% Chlorhexidine CAS 55-56-1  
≤10% Didecyl Dimethyl Ammonium chloride CAS 7173-51-5  
≤0.5% Quaternary Ammonium Compounds CAS 68956-79-6  
<5% Benzyl-C 12-18 alkyldimethylammonium chlorides CAS 63891-01-5

#### Others:

≤10% Glutaral CAS 111-30-8  
≤0.63% Sodium Hypochlorite CAS 7681-52-9  
≤7.5% Povidone Iodine CAS 25655-41-8  
0.5% Disodium Cocoampho Dipropionate CAS68604-71-7  
≤1.4% Hydrogen Peroxide CAS 7722-84-1

Access the latest Approved Cleaners list via [adeptmedicaltraining.com/downloads](http://adeptmedicaltraining.com/downloads)

## Disposal

Used product is a biohazard, decontaminate according to instruction provided in the Disinfection Instructions and reach out to the manufacturer on [adeptmedical@adept.co.nz](mailto:adeptmedical@adept.co.nz) for material specification, if required. Dispose according to internal clinic policy taking into consideration local regulations.

[adeptmedical@adept.co.nz](mailto:adeptmedical@adept.co.nz)

### Maintenance, Inspection and Testing

#### For all components ensure:

- No visible damage or sharp edges, e.g. cracks.

#### Retrograde Main Platform, Front Base, Extension (when applicable):

- No exposed Carbon Fibre is visible.

#### For the Extension (when applicable):

- Ensure the Extension can be assembled smoothly onto the Retrograde Main Platform and the locking clips are functional.

#### For the Front Base and Rear Base:

- Ensure the locking clips are functional and can retain the Leg.

#### For the Leg:

- Ensure thread is intact and can be screwed onto the Retrograde Main Platform.
- Ensure the blue collar rotates smoothly.
- Ensure when blue collar is in the unlocked position the Leg can extend and contract smoothly when force is applied.
- Ensure when blue collar is in the locked position the Leg cannot extend and contract when force is applied.

#### For the Leg Extension:

- Ensure thread is intact and can be screwed onto the Leg Extension and the Retrograde Main Platform.

### Packaging

It is not required to package the device following disinfection.

### Sterilisation

This device should not be subjected to sterilisation processes.

### Storage

Once disinfection is complete and all components are dry the device should be stored in a dry environment.

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

### Serious Incidents

Any serious incident which occurs in relation to the device should be reported immediately to:

- The Competent Authority of the applicable Member State.
- The Manufacturer using the Serious Incident Reporting Form found on the Adept Medical Website.

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